

PAX4 Rabbit pAb

货号: B17047

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa/37kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Lovo
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human PAX4 (NP_006184.2).
序列	LGRYYRTGVLEPKGIGGSKPRLATPPVVARIAQLKGECPALFAWEIQRQLCAEGLCTQDKTPSVSSINRVLRALQEDQGLPC TRLRSPAVLAPAVLTPHSGSETPRGTHPGTGHRNRTIFSPSQAEALEKEFQRGQYPDSVARGKLATATS

靶点信息

研究背景	This gene is a member of the paired box (PAX) family of transcription factors. Members of this gene famil y typically contain a paired box domain, an octapeptide, and a paired-type homeodomain. These genes pl ay critical roles during fetal development and cancer growth. The paired box 4 gene is involved in pancre atic islet development and mouse studies have demonstrated a role for this gene in differentiation of insu lin-producing beta cells.
基因ID	5078

基因名	PAX4
Swiss	O43316
别名	PAX4;KPD;MODY9

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn