

DHCR24 Rabbit pAb

货号: B17055

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	54kDa/60kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HT-29,PC-12,NCI-H460
细胞定位	Endoplasmic reticulum membrane,Golgi apparatus membrane,Single-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 227-516 of human DHC R24 (NP_055577.1).
序列	AEIRIIPAKKYVKLRFEPVRGLEAICAKFTHESQRQENHFVEGLLYSLDEAVIMTGVMTDEAEPKLNLSIGNYYKPWFVKHVE NYLKTNREGLEYIPLRHYYHRHTRSIFWELQDIIPFGNNPIFRYLFGWMPKISLLKLTQGETLRKLYEQHHVVQDMLVPMK CLQQALHTFQNDIHVYPIWLCPFILPSQPGLVHPKGNEAELYIDIGAYGEPRVKHFEARSCMRQLEKFVRSVHGFQMLYAD CYMNREEFWEMFDGSLYHKLREKLGQCDAFPEVYDKICKAARH

靶点信息

研究背景	This gene encodes a flavin adenine dinucleotide (FAD)-dependent oxidoreductase which catalyzes the red uction of the delta-24 double bond of sterol intermediates during cholesterol biosynthesis. The protein co ntains a leader sequence that directs it to the endoplasmic reticulum membrane. Missense mutations in t his gene have been associated with desmosterolosis. Also, reduced expression of the gene occurs in the t emporal cortex of Alzheimer disease patients and overexpression has been observed in adrenal gland ca ncer cells.
------	---

基因ID	1718
基因名	DHCR24
Swiss	Q15392
别名	DHCR24;DCE;Nbla03646;SELADIN1;seladin-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn