

Alpha2-adaptin/AP2A2 Rabbit pAb

货号: B17060

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	73kDa/103kDa/104kDa
实测分子量	104KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse liver,Rat liver
细胞定位	Cell membrane,Cytoplasmic side,Membrane,Peripheral membrane protein,coated pit
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 590-939 of human Alph a2-adaptin/AP2A2 (NP_036437.1).
序列	STDILATVLEEMPPFPERESSILAKLKKKKGPSTVTDLEDTKRDRSVDVNGGPEPAPASTSAVSTPSPSADLLGLGAAPPAP AGPPPSSGGSGLLVDVFSDSASVVAPLAPGSEDNFARFVCKNNGVLFENQLLQIGLKSEFRQNLGRMFIFYGNKTSTQFL NPTPTLICSDDLQPNLNLQTKPVDPTVEGGAQVQQVWNIECVSDFTEAPVLNIQFRYGGTFQNVSVQLPITLNKFFQPTM ASQDFFQRWKQLSNPQQEVQNIFKAKHPMDTEVTAKIIGFGSALLEEVDPNPANFVGAGIIHTKTTQIGCLLRLEPNLQA QMYRLTLRTSKEAVSQRLCELLSAQF

靶点信息

研究背景	The protein encoded by this gene is a subunit of the AP-2 adaptor protein complex, which is involved in li nking lipid and protein membrane components with the clathrin lattice. This interaction supports the form ation of clathrin-coated vesicles, and the encoded subunit aids in the process by binding polyphosphoinos itide-containing lipids in the cell membrane.
------	--

基因ID	161
基因名	AP2A2
Swiss	O94973
别名	AP2A2;ADTAB;CLAPA2;HIP-9;HIP9;HYPJ

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn