

NEU4 Rabbit pAb

货号: B17112

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa/52kDa/53kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,Mouse pancreas,Mouse liver
细胞定位	Lysosome lumen,Membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 197-496 of human NEU 4 (NP_542779.2).
序列	ECFGKICRTSPHSFAFYSDDHGRTWRCGGLVPNLRSGECLAAVDGGQAGSFLYCNARSPLGSRVQALSTDEGTSFLPA ERVASLPETAWGCQGSIVGFAPAPNRPDDSWSVGPGSPLQPPLLGPGVHEPPEEAAVDPRGGQVPGGPF SRLQPRG DGPRQGP RGPVSGDVGSWTLALPMPFAAPPQSPTWLLYSHPVGRRARLHMGIRLSQSPLDPRSWTEPWVIYEGPSGY SDLASIGPAPEGGLVFACLYESGARTSYDEISFCTFSLRELVNPASPKPPNLGDKPRGCCWPS

靶点信息

研究背景	The protein encoded by this gene belongs to a family of glycohydrolytic enzymes, which remove terminal sialic acid residues from various sialo derivatives, such as glycoproteins, glycolipids, oligosaccharides, and gangliosides. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
------	--

基因ID	129807
基因名	NEU4
Swiss	Q8WWR8
别名	NEU4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn