

LARP4 Rabbit pAb

货号: B17114

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	49kDa/67kDa/72kDa/80kDa/81kDa
实测分子量	81KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat kidney
细胞定位	cytoplasmic stress granule,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human LARP4 (NP_443111.4).
序列	MLLFVEQVASKGTGLNPNAKVWQEIAPGNTDATPVTHGTESSWHEIAATSGAHPEGNAELSEDICKEYEVMYSSSCETTR NTTGIEESTDGMILGPEDLSYQIYDVSGESNSAVSTEDLK

靶点信息

研究背景	RNA binding protein that binds to the poly-A tract of mRNA molecules. Associates with the 40S ribosomal subunit and with polysomes. Plays a role in the regulation of mRNA translation. Plays a role in the regulation of cell morphology and cytoskeletal organization.
基因ID	113251

基因名	LARP4
Swiss	Q71RC2
别名	LARP4;PP13296

产品验证



Western blot analysis of LARP4 expressed in Rat kidney using LARP4 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn