

ATRIP Rabbit pAb

货号: B17116

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	72kDa/76kDa/83kDa/85kDa
实测分子量	100kDa/113kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ATRIP (NP_569055.1).
序列	MAGTSAPGSKRRSEPPAPRPGPPPGTGHPPSKRARGFSAAAAPDPDDPFGAHGDFTADDLEELDTLASQALSQCPAAAR DVSSDHKVRLLDGM SKNPSGKNRETVP IKDNFELEVLQAQYKELKEKMKVMEEEVLIKNGEIKILRDSLHQTESVLEEQR RSHFLLEQEK TQALSDKEKEFSKKLQSLQSELQFKDAEMNELRTKLQTSERANKLAAPSVSHVSPRKNPSVVIKPEACSPQ FGKTSFPTKESFSANMSLP

靶点信息

研究背景	This gene encodes an essential component of the DNA damage checkpoint. The encoded protein binds to single-stranded DNA coated with replication protein A. The protein also interacts with the ataxia telangieci tasia and Rad3 related protein kinase, resulting in its accumulation at intranuclear foci induced by DNA d amage. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	84126
基因名	ATRIP
Swiss	Q8WXE1
别名	ATRIP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn