

CHMP6 Rabbit pAb

货号: B17121

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	23kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U2OS,PC-3,Jurkat
细胞定位	Endomembrane system,Endosome membrane,Late endosome membrane,Lipid-anchor,Membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 111-200 of human CHMP6 (NP_078867.2).
序列	LNKMHQVMSIEEVERILDETQEAHEYQRQIDELLAGSFTQEDEDAILLELSAITQEIQIELPEVPSEPLPEKIPENVPVKARPRQAELVAA

靶点信息

研究背景	This gene encodes a member of the chromatin-modifying protein/charged multivesicular body protein family. Proteins in this family are part of the ESCRT-III (endosomal sorting complex required for transport III) which degrades surface receptors, and in biosynthesis of endosomes.
基因ID	79643

基因名	CHMP6
Swiss	Q96FZ7
别名	CHMP6;VPS20

产品验证



Western blot analysis of CHMP6 expressed in U2OS,PC-3,Jurkat using CHMP6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn