

WDFY4 Rabbit pAb

货号: B17129

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:100
理论分子量	73kDa/114kDa/325kDa/352kDa/353kDa
实测分子量	354kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human WDFY4 (NP_065996.1).
序列	MEAEDLSKAEDRNEDPGSKNEGQLAAVQPDVPHGGQSSSPTALWDMLEKRFLEYQQLTHKSPIERQKSLLSLLPLFLKA WEHSVGIICFSLQRLAEDVSDQLAQLQKALVGKPAEQARLAAGQLLWWKGDVDQDGYLLLSVYVLTGTDSETLGRV AE

靶点信息

研究背景	Plays a critical role in the regulation of cDC1-mediated cross-presentation of viral and tumor antigens in d endritic cells. Mechanistically, acts near the plasma membrane and interacts with endosomal membranes to promote endosomal-to-cytosol antigen trafficking. Plays also a role in B-cell survival through regulation of autophagy.
------	---

基因ID	57705
基因名	WDFY4
Swiss	Q6ZS81
别名	WDFY4;C10orf64

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn