

RHOF Rabbit pAb

货号: B17141

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	18kDa/23kDa
实测分子量	27kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,293T,Mouse kidney
细胞定位	Cell membrane,Cytoplasm,Cytoplasmic side,Lipid-anchor,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human RHOF (NP_061907.2).
序列	SQGSFPEHYAPSVFEKYTASVTVGSKVETLNLTYDTAGQEDYDRLRPLSYQNTHLVLCYDVMNPTSNDVNLIKWFPVTHF

靶点信息

研究背景	Plasma membrane-associated small GTPase which cycles between an active GTP-bound and an inactive G DP-bound state. Causes the formation of thin, actin-rich surface projections called filopodia. Functions co operatively with CDC42 and Rac to generate additional structures, increasing the diversity of actin-based morphology.
基因ID	54509

基因名	RHOF
Swiss	Q9HBH0
别名	RHOF;ARHF;RIF

产品验证



Western blot analysis of RHOF expressed in HT-29,293T,Mouse kidney using RHOF Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Block ing buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn