

MYCBP Rabbit pAb

货号: B17155

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	11kDa
实测分子量	12kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,K-562
细胞定位	Cytoplasm,Mitochondrion,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human MYCBP (NP_036465.2).
序列	MAHYKAADSKREQFRRYLEKSGVLDLTKVLVALYEEPEKPNsALDFLKHHLGAATPENPEIELLRLELAEMKEKYEAIVEEN KKLKAKLAQYEPPQEEKRAE

靶点信息

研究背景	The protein encoded by this gene binds to the N-terminus of the oncogenic protein C-MYC, enhancing the ability of C-MYC to activate E box-dependent transcription. The encoded protein is normally found in the cytoplasm, but it translocates to the nucleus during S phase of the cell cycle and associates with C-MYC. This protein may be involved in spermatogenesis. This gene can be silenced by microRNA-22. Two transcript variants, one protein-coding and the other probably not protein-coding, have been found for this gene.
------	--

基因ID	26292
基因名	MYCBP
Swiss	Q99417
别名	MYCBP;AMY-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn