

# BLOC1S6 Rabbit pAb

货号: B17156

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB                                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:200 - 1:2000                                                                              |
| 理论分子量 | 8kDa/19kDa                                                                                             |
| 实测分子量 | 20kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | U-87MG,HepG2,DU145,K-562,Rat testis                                                                    |
| 细胞定位  | Cytoplasm,Endomembrane system,Peripheral membrane protein                                              |
| 纯化    | Affinity purification                                                                                  |

抗原信息


|      |                                                                                                                                                                            |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human BLOC1 S6 (NP_036520.1).                                                       |
| 序列   | MSVPGPSSPDGALTRPPYCLEAGEPTPGLSDTSPDEGLIEDLTIEDKAVEQLAEGLLSHYLPDLQRSKQALQELTQNQVVLLDTLEQEISKFKECHSMLDINALFAEAKHYHAKLVNIRKEMLMLEKTSKLLKKRALKLQQKRQKEELEREQQREKEFEREKQLTARPAK |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                         |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene may play a role in intracellular vesicle trafficking. It interacts with Synta xin 13 which mediates intracellular membrane fusion. Mutations in this gene cause symptoms associated with Hermansky-Pudlak syndrome-9. Alternative splicing results in multiple transcript variants. A pseudog ene related to this gene is located on the X chromosome. |
| 基因ID | 26258                                                                                                                                                                                                                                                                                                                                                                                   |

|       |                                   |
|-------|-----------------------------------|
| 基因名   | BLOC1S6                           |
| Swiss | Q9UL45                            |
| 别名    | BLOC1S6;BLOS6;HPS9;PA;PALLID;PLDN |

产品验证



Western blot analysis of BLOC1S6 expressed in Rat testis,U-87MG,HepG2 using BLOC1S6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)