

PPIL2 Rabbit pAb

货号: B17162

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	58kDa/59kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse spleen
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PPIL2 (NP_680481.1).
序列	MGKRQHQKDKMYITCAEYTHFYGGKKPDLPQTNFRRLPFDHCSSLQPFVYPVCTPDGIVFDLLNIVPWLKKYGTNPNSNG EKLDGRSLIKLNFSKNSEGKYHCPVLFTVFTNNTNTHIVAVRTTGNVYAYEAVEQLNIKAKNFRDLLTDEPFSRQDIITLQDPTN LDKFNVSNFYHVKNMKIIDPDEEKAKQDPSYYLKNTNAETRETQLQELYKEFKGDEILAATMKAPEKKKVDKLNAAHYSTG KVSASFTSTAMPETT

靶点信息

研究背景	This gene is a member of the cyclophilin family of peptidylprolyl isomerases. The cyclophilins are a highly conserved ubiquitous family, members of which play an important role in protein folding, immunosuppression by cyclosporin A, and infection of HIV-1 virions. This protein interacts with the proteinase inhibitor eg lin c and is localized in the nucleus. Alternatively spliced transcript variants have been found for this gene .
------	---

基因ID	23759
基因名	PPIL2
Swiss	Q13356
别名	PPIL2;CYC4;CYP60;Cyp-60;UBOX7;hCyP-60

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn