

ADNP Rabbit pAb

货号: B17170

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	123kDa
实测分子量	163kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-87MG,HT-29,Mouse liver
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 843-1102 of human ADNP (NP_056154.1).
序列	EWLFENHDEKDSRVNASKTADKKLNLGKEDDSSSDSFENLEEEESNESGSPFDPVFEVEPKISNDNP EEHVLKVIPEDASE SEEKLDQKEDGSKYETIHLTEEPTKLMHNASDSEVDQDDVVEWKDGASPS ESGPGSQVSDFEDNTCEMKPGTWSDE SSQSEDARSSKPAAKKKATMQGDREQLKWKNSSYGKVEGFWSKDQSQWKNASENDERLSNPQIEWQNSTIDSEGE QFDNMTDGVAEPMHGSLAGVKLSSQQA

靶点信息

研究背景	Vasoactive intestinal peptide is a neuroprotective factor that has a stimulatory effect on the growth of some tumor cells and an inhibitory effect on others. This gene encodes a protein that is upregulated by vasoactive intestinal peptide and may be involved in its stimulatory effect on certain tumor cells. The encoded protein contains one homeobox and nine zinc finger domains, suggesting that it functions as a transcription factor. This gene is also upregulated in normal proliferative tissues. Finally, the encoded protein may increase the viability of certain cell types through modulation of p53 activity. Alternatively spliced transcript variants encoding the same protein have been described.
基因ID	23394
基因名	ADNP
Swiss	Q9H2P0
别名	ADNP;ADNP1;HVDAS;MRD28

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn