

CRTAP Rabbit pAb

货号: B17194

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	46kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human CRTAP (NP_006362.1).
序列	ALRAGRAQYERYSFRRDELMPLESAYRHALDKYSGEHWAESVGYLEISLRLHRLLRDSEAFCHRNCSAAPQPEPAAGL ASYPELRLFGGLLRRRAHCLKRCKQGLPAFRQSQPSREVLADFQRREPYKFLQFAYFKANNLPKAIAAAHTFLLKHPDDEMM KRNMAYYKSLPGAEDYIKDLETKSYESLFIRAVRAYNGENWRTSITDMELALPDFFKAF

靶点信息

研究背景	The protein encoded by this gene is similar to the chicken and mouse CRTAP genes. The encoded protein is a scaffolding protein that may influence the activity of at least one member of the cytohesin/ARNO family in response to specific cellular stimuli. Defects in this gene are associated with osteogenesis imperfecta, a connective tissue disorder characterized by bone fragility and low bone mass.
------	--

基因ID	10491
基因名	CRTAP
Swiss	O75718
别名	CRTAP;CASP;LEPREL3;OI7;P3H5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn