

EDIL3 Rabbit pAb

货号: B17199

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	52kDa/53kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,DU145
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-300 of human EDIL3 (NP_005702.3).
序列	SSVVEVASDEEEPTSGAPCTPNPCHNGGTCEISEAYRGDTFIGYVCKCPRGFNGIHCQHNECEVEPCKNGGICTDLVAN YSCECPGEFMGRNCQYKCSGPLGIEGGIISNQITASSTHRALFGLQKWYPYARLNKKGLINAWTAAENDRWPWIQINLQ RKMRVTGVITQGAKRIGSPEYIKSYKIAYSNDGKTWAMYKVKGTTNEDMVFRGNIDNNTPYANSFTPPIKAQYVRLYPQV

靶点信息

研究背景	The protein encoded by this gene is an integrin ligand. It plays an important role in mediating angiogenes is and may be important in vessel wall remodeling and development. It also influences endothelial cell be havior.
基因ID	10085

基因名	EDIL3
Swiss	O43854
别名	EDIL3;DEL1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn