

# UXT Rabbit pAb

货号: B17220

## 产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                  |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB                                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                              |
| 理论分子量 | 18kDa                                                                                                  |
| 实测分子量 | 17kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | 293T,HL-60                                                                                             |
| 细胞定位  | Cytoplasm,Nucleus,centrosome,cytoskeleton,microtubule organizing center                                |
| 纯化    | Affinity purification                                                                                  |

## 抗原信息


|      |                                                                                                                                                                 |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 13-169 of human UXT ( NP_705582.1).                                               |
| 序列   | MATPPKRRAVEATGEKVLRYETFISDVLQRDLRKVLDHRDKVYEQLAKYLQLRNVIERLQEAKHSELYMQVDLGCNFFVDT<br>VVPDTSRIYVALGYGFFLELTAEALKFIDRKSSLLTELSNSLTKD SMNIAHIMLLEGLRELQGLQNFPEKPHH |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent tr anscription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this ge ne may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms h ave been observed for this gene. |
| 基因ID | 8409                                                                                                                                                                                                                                                                                                                                                              |

|       |                  |
|-------|------------------|
| 基因名   | UXT              |
| Swiss | Q9UBK9           |
| 别名    | UXT;ART-27;STAP1 |

产品验证



Western blot analysis of UXT expressed in 293T,HL-60 using UXT Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)