

RBM10 Rabbit pAb

货号: B17221

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	94kDa/103kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,MCF7,LO2,293T
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human RBM10 (NP_005667.2).
序列	MEYERRGGRGDRTGRYGATDRSQDDGGENRSRDHDYRDMDYRSYPREYGSQEGKHDYDDSSSEEQSAEDSYEASPGSETQRRRRRRHRHSPTGPPGFPRDGDYRDQDYRTEQEEEEEEEEEEEEKASNIVMLRMLPQAATEDDIRGQLQSHGVQAREVRLMRNKSSGQSRGFAFVEFSH

靶点信息

研究背景	This gene encodes a nuclear protein that belongs to a family proteins that contain an RNA-binding motif. The encoded protein associates with hnRNP proteins and may be involved in regulating alternative splicing. Defects in this gene are the cause of the X-linked recessive disorder, TARP syndrome. Alternate splicing results in multiple transcript variants.
基因ID	8241

基因名	RBM10
Swiss	P98175
别名	RBM10;DXS8237E;GPATC9;GPATCH9;S1-1;TARPS;ZRANB5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn