

FZD7 Rabbit pAb

货号: B17223

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	63kDa
实测分子量	70-72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,Mouse kidney
细胞定位	Cell membrane,Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human FZD7 (NP_003498.1).
序列	MRDPGAAAPLSSLGLCALVLALLGALSAGAGAQPYHGEKGISVPDHGFCQPISIPLCTDIAYNQITLPNLLGHTNQEDAGLEVHQFYPLVKVQCSPFLRFFLCSEMYAPVCTVLDQAIPPCRSQCERARQGCEALMNKFGFWPERLRNCFPVHGAGEICVGGQNTSDGSGGPGGGPTAYPTAPYLPDLPTALPPGASDGRGRPAFPFSCPRQLKVPYLYGYRFLGERDCGAPCEPGRANGLMYFKEEERRFARLWVGWV

靶点信息

研究背景	Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD7 protein contains an N-terminal signal sequence, 10 cysteine residues typical of the cysteine-rich extracellular domain of Fz family members, 7 putative transmembrane domains, and an intracellular C-terminal tail with a PDZ domain-binding motif. FZD7 gene expression may downregulate APC function and enhance beta-catenin-mediated signals in poorly differentiated human esophageal carcinomas.
------	--

基因ID	8324
基因名	FZD7
Swiss	O75084
别名	FZD7;FzE3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn