

ARHGEF5 Rabbit pAb

货号: B17225

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:200
理论分子量	60kDa/176kDa
实测分子量	176kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,A-549,HepG2,293T,Rat lung
细胞定位	cell projection,cytoplasm,cytosol,nucleoplasm,nucleus,plasma membrane,podosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1248-1597 of human ARHGEF5 (NP_005426.2).
序列	QVCDVVLNHAPDFRRVYLPYVTNQTYQERTFQSLMNSNSNFREVLEKLESDPVCQRLSLKSFLILPFQRITRLKLLLQNILKR TQPGSSEAEATKAHHALEQLIRDCNNNVQSMRRTTEELIYLSQKIEFECKIFPLISQSRWLKSGELTALEFSASPGLRRKLN TRPVHLHLFNDCLLSRPREGSRFLVFDHAPFSSIRGEKCEMKLHGPHKNLFRFLRQNTQGAQAEFLFRTETQSEKLRWIS ALAMPREELDLLECYNPQVQCLRAYKPRENDELALEKADVVMVTQQSSDGWLEGVRLSDGERGWFPVQQVEFISNPEV RAQNLKEAHRVKTAQLLVEQQA

靶点信息

研究背景	Rho GTPases play a fundamental role in numerous cellular processes initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein may be involved in the control of cytoskeletal organization.
------	--

基因ID	7984
基因名	ARHGEF5
Swiss	Q12774
别名	ARHGEF5;GEF5;P60;TIM;TIM1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn