

TMOD1 Rabbit pAb

货号: B17232

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:4000
理论分子量	26kDa/40kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,Mouse heart,Mouse skeletal muscle,Mouse lung
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-359 of human TMOD1 (NP_003266.1).
序列	MSYRRELEKYRDLDEDEILGALTEELRTLENELDELDPDNALLPAGLRQKDQTTKAPTGPFKREELLDHLEKQAKEFKDRE DLVPYTGEKRGKVWVPKQKPLDPVLESVTLEPELEEALANASDAELCDIAAILGMHTLMSNQYYQALSSSSIMNKEGLNS VIKPTQYKVPDEEPNSTDVEETLERIKNNDPKLEEVNLLNIRNIPIPTLKAYAEALKENSYVKKFSIVGTRSNDPVAYALAEML KENKVLKTLNVESNFISGAGILRLVEALPYNTSLVEMKIDNQSQPLGNKVEMEIVSMLEKNATLLKFGYHFTQQGPRLRASN AMMNNNDLVRKRRRLADLTGPIPKCRSGV

靶点信息

研究背景	This gene encodes a member of the tropomodulin family. The encoded protein is an actin-capping protein that regulates tropomyosin by binding to its N-terminus, inhibiting depolymerization and elongation of the pointed end of actin filaments and thereby influencing the structure of the erythrocyte membrane skeleton. Multiple transcript variants encoding the same protein have been found for this gene.
------	--

基因ID	7111
基因名	TMOD1
Swiss	P28289
别名	TMOD1;D9S57E;ETMOD;TMOD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn