

SH3GL3 Rabbit pAb

货号: B17238

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:100 - 1:200
理论分子量	31kDa/38kDa/39kDa/40kDa
实测分子量	39KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cytoplasm,Early endosome membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human SH3GL3 (NP_003018.3).
序列	MSVAGLKKQFHKASQLFSEKISGAEGTKLDDEFLDMERKIDVTNKVVAEILSKTTEYLQPNPAYRAKLGMLNTVSKIRGQV KTTGYPQTEGLLGDCMLKYGKELGEDSTFGNALIEVGESMKLMAEVKDSLINVKQTFIDPLQLLQDKDLKEIGHHLKKLE GRRLDYDYKKKRVGKIPDEEVRQAVEKFEEKELAERSMFNPLENDVEQVSQLAVFIEAALDYHRQSTEILQELQSKLQMRI SAASSVPRREYKPRPVKRSSSELNGV

靶点信息

研究背景	Implicated in endocytosis. May recruit other proteins to membranes with high curvature (By similarity.
基因ID	6457

基因名	SH3GL3
Swiss	Q99963
别名	SH3GL3;CNSA3;EEN-B2;HsT19371;SH3D2C;SH3P13

产品验证



Western blot analysis of SH3GL3 expressed in Mouse brain,Rat brain using SH3GL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn