

QSOX1 Rabbit pAb

货号: B17243

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	66kDa/82kDa
实测分子量	66KFa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460,U-251MG,Mouse heart
细胞定位	Golgi apparatus membrane,Secreted,Single-pass membrane protein,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human QSOX1 (NP_001004128.1).
序列	MRRCNSGSGPPPSLLLLLLWLLAVPGANAAPRSALYSPSDPLTLQADTVRGAVLGSRSAWAVEFFASWCGHCIAFAPTW KALAEDVKAWRPALYLAALDCAEETNSAVCRDFNIPGFPTVRFFKAFTKNGSGAVFPVAGADVQTLRERLIDALESHHDT WPPACPPLEPAKLEEIDGFFARNNEEYLALIFEKGGSYLGREVALDLSQHKGVAVRRVLNTEANVVRKFGVTDFPSCYLLFR NGSVSRVPVLMESRSFYTAYLQRLSGLTREAAQTTVAPTANKIAPT VWKLADRSKIYMADLESALHYILRIEVGRFPVLEGQ RLVAL

靶点信息

研究背景	This gene encodes a protein that contains domains of thioredoxin and ERV1, members of two long-standing gene families. The gene expression is induced as fibroblasts begin to exit the proliferative cycle and enter quiescence, suggesting that this gene plays an important role in growth regulation. Two transcript variants encoding two different isoforms have been found for this gene.
------	---

基因ID	5768
基因名	QSOX1
Swiss	O00391
别名	QSOX1;Q6;QSCN6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn