

GPR17 Rabbit pAb

货号: B17257

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	37kDa/40kDa
实测分子量	41KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spinal cord,Mouse brain,Rat brain
细胞定位	plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-94 of human GPR17 (NP_005282.1).
序列	MSKRSWWAGSRKPPREMLKLSGSDSSQSMNGLEVAPPGLITNFSLATAEQCGQETPLENMLFASFYLLDFILALVGNTLA LWLFIRDHKSGTPA

靶点信息

研究背景	Dual specificity receptor for uracil nucleotides and cysteinyl leukotrienes (CysLTs. Signals through G(i and inhibition of adenylyl cyclase. May mediate brain damage by nucleotides and CysLTs following ischemia.
基因ID	2840
基因名	GPR17

Swiss	Q13304
别名	GPR17

产品验证



Western blot analysis of GPR17 expressed in Mouse spinal cord,Mouse brain,Rat brain using GPR17 Rab bit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30 ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn