

GRK5 Rabbit pAb

货号: B17258

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	68kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,U-251MG,HepG2,293T,HeLa,HT-1080,Mouse liver
细胞定位	Cell membrane,Cytoplasm,Nucleus,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 371-590 of human GRK 5 (NP_005299.1).
序列	LGCLYEMIEGQSPFRGRKEKVKREEVDRRVLETEEVYSHKFSEEAKSICKMLLTKDAKQRLGCQEEGAAEVKRHPFFRNM NFKRLEAGMLDPPFVPDPRAVYCKDVLDIEQFSTVKGVNLDHTDDDFYSKFSTGVSIPWQNEMIETECFKELNVFGPNG TLPPDLNRNHPPEPPKGLLQRLFKRQHQNNSKSSPSSKTSFNHHINSNHVSSNSTGSS

靶点信息

研究背景	This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled receptor kinase subfamily of the Ser/Thr protein kinase family. The protein phosphorylates the activated forms of G protein-coupled receptors thus initiating their deactivation. It has also been shown to play a role in regulating the motility of polymorphonuclear leukocytes (PMNs).
------	---

基因ID	2869
基因名	GRK5
Swiss	P34947
别名	GRK5;GPRK5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn