

CAPZA1 Rabbit pAb

货号: B17280

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:1000
理论分子量	32kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,MCF7,HL-60,Mouse brain,Mouse kidney,Mouse lung,Rat lung
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-286 of human CAPZA1 (NP_006126.1).
序列	MADFDDRVSDEEKVRIAAKFITHAPPGEFNEVFNDVRLLLNNDNLLREGAAHAFAQYNMDQFTPVKIEGYEDQVLITEHG DLGNSRFLDPRNKISFKFDHLRKEASDPQPEEADGGLKSWRESCDSALRAYVKDHYSNGFCTVYAKTIDGQQTIACIESH QFQPKNFWNGRWRSEWKFTITPPTAQVVGVLKIQVHYIEDGNVQLVSHKDVQDSLTVSNEAQTAKFIKIENAENEYQT AISENYQTMSDITTFKALRRQLPVTRTKIDWNKILSYKIGKEMQNA

靶点信息

研究背景	CAPZA1 is a member of the F-actin capping protein alpha subunit family. This gene encodes the alpha subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments.
基因ID	829

基因名	CAPZA1
Swiss	P52907
别名	CAPZA1;CAPPA1;CAPZ;CAZ1

产品验证



Western blot analysis of CAPZA1 expressed in Mouse brain,Rat lung,HeLa using CAPZA1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn