

TRIM41 Rabbit pAb

货号: B17303

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	23kDa/43kDa/59kDa/71kDa
实测分子量	72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,LO2,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 421-630 of human TRIM 41 (NP_291027.3).
序列	MSRMFCQAARVDLTLPDTAHPALMLSPDRRGVRLAERRQEVAADHPKRFSADCCVLGAQGFRSGRHYWEVEVGRRG WAVGAARESTHHKEKVGPGGSSVSGDASSRHHHRRRLHLPQQPLLQREVWCVTNGKRYQAQSSTEQTLLSPSE KPRRFGVYLDYEAGRLGFYNAETLAHVHTFSAAFLGERVFPFFRVL SKGTRIKLCP

靶点信息

研究背景	This gene encodes a member of the tripartite motif (TRIM) family. The TRIM family is characterized by a signature motif composed of a RING finger, one or more B-box domains, and a coiled-coil region. This encoded protein may play a role in protein kinase C signaling. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	90933
基因名	TRIM41
Swiss	Q8WV44
别名	TRIM41;RINCK

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn