

AFAP1 Rabbit pAb

货号: B17309

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:1000 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	80kDa/89kDa
实测分子量	81kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,A-549,HeLa,U-251MG,HepG2,Mouse liver,Rat liver
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 461-730 of human AFAP 1 (NP_940997.1).
序列	SVIQTAKQTF CFMNRVISANPYLGGTSNGYAHPSGTALHYDDVPCINGSLKGKKPPVASNGVTGKGKTLSSQPKKADPA AVVKRTGSNAAQYKYGKNRVEADAKRLQTKEEELLKRKEALRNRLAQLRKERKDLRAAIEVNAGRKPQAILEEKLKQLEEEC RQKEAERSVLELEL TEVKESLKKALAGGVTLGLAIEPKSGTSSPQSPVFRHRTLENSPISSCDTSDTEGPVPVNSAAVLKKS QAAPGSSPCRGHVLRKAKEWELKNGT

靶点信息

研究背景	The protein encoded by this gene is a Src binding partner. It may represent a potential modulator of actin filament integrity in response to cellular signals, and may function as an adaptor protein by linking Src family members and/or other signaling proteins to actin filaments. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	60312
基因名	AFAP1
Swiss	Q8N556
别名	AFAP1;AFAP;AFAP-110;AFAP110

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn