

CREBZF Rabbit pAb

货号: B17310

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	37kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,293T,HeLa,U-87MG,K-562,Mouse brain,Mouse heart,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-354 of human CREBZF (NP_001034707.1).
序列	MRHSLTKLLAASGSNSPTRSESPAPAATCSLPDLTRAAAGEEETAAAGSPGRKQQFGDEGELEAGRGSRGGVAVRAPSP PEEMEEEAIASLPGEETEDMDFLSGLELADLLDPRQPDWHLDPGLSSPGPLSSSGGGSDSGGLWRGDDDDDEAAAAEMQ RFSDLLQRLNLNGIGGCSSSSDSGSAEKRRRKSPGGGGGGSGNDNNQAATKSPRKAATAAARLNRLKKKEYVMGLES VRGLAAENQELRAENRELGKRVQALQEESRYLRANLANETGLARLLSRLSGVGLRLTTSLFRDSPAGDHDYALPVGKQKQ DLLEEDDSAGGVCLHVDKDKVSVEFCSACARKASSSLKM

靶点信息

研究背景	Strongly activates transcription when bound to HCFC1. Suppresses the expression of HSV proteins in cells infected with the virus in a HCFC1-dependent manner. Also suppresses the HCFC1-dependent transcriptional activation by CREB3 and reduces the amount of CREB3 in the cell. Able to down-regulate expression of some cellular genes in CREBZF-expressing cells.
------	--

基因ID	58487
基因名	CREBZF
Swiss	Q9NS37
别名	CREBZF;SMILE;ZF

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn