

COX8A Rabbit pAb

货号: B17337

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	7kDa
实测分子量	8kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse skeletal muscle
细胞定位	Mitochondrion inner membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-69 of human COX8A (NP_004065.1).
序列	MSVLTPLLLRGLTGSARRLPVPRAKIHSLPPEGKLGIMELAVGLTSCFVTFLLPAGWILSHLETYRRPE

靶点信息

研究背景	The protein encoded by this gene is the terminal enzyme of the respiratory chain, coupling the transfer of electrons from cytochrome c to molecular oxygen, with the concomitant production of a proton electroch emical gradient across the inner mitochondrial membrane. In addition to 3 mitochondrially encoded subu nits, which perform the catalytic function, the eukaryotic enzyme contains nuclear-encoded smaller subu nits, ranging in number from 4 in some organisms to 10 in mammals. It has been proposed that nuclear-e ncoded subunits may be involved in the modulation of the catalytic function. This gene encodes one of th e nuclear-encoded subunits.
------	--

基因ID	1351
基因名	COX8A
Swiss	P10176
别名	COX8A;COX;COX8;COX8-2;COX8L;VIII;VIII-L

产品验证



Western blot analysis of COX8A expressed in Mouse skeletal muscle using COX8A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn