

NFIA Rabbit pAb

货号: B17341

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	54kDa/55kDa/60kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 180-330 of human NFIA (NP_005586.1).
序列	AYFVHAADSSQSESPSQPSDADIKDPENHGLGFQDSFVTSGVFSVTELVVRVSQTPIAAGTGPNFSLSDLESSSYYSMSP GAMRRSLPSTSTSTKRLKSVEDEMDSPGEEPFTGTQGRSPGSGSQSSGWHEVEPGMPSPPTTLKKSEKSG

靶点信息

研究背景	This gene encodes a member of the NF1 (nuclear factor 1) family of transcription factors. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	4774

基因名	NFIA
Swiss	Q12857
别名	NFIA;CTF;NF-I/A;NF1-A;NFI-A;NFI-L;BRMUTD

产品验证



Western blot analysis of NFIA expressed in Mouse brain using NFIA Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn