

NELFE Rabbit pAb

货号: B17344

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	40kDa/43kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,HepG2,SW480,Mouse liver,Mouse testis,Mouse thymus
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NELFE (NP_002895.3).
序列	MLVIPPLSEEEEEALQKKFNKLKKKKKALLALKKQSSSSTTSQGGVKRSLSEQPVMDTATATEQAKQLVKSGAISAIAETK NSGFKRSRTLEGKLDPEKGPVPTFQPFQRSISADDDLQESSRRPQRKSLYESFVSSSDRLRELGPDGEEAEGPGAGDGP PRSFWDWGYEERSGAHSSASPPRSRSRDRSHERNRDRDRDRERDRDRDRDRDRERDRDRDRDRDRERDRDRERDR DRDREGPFRRSDSFPERRAPRKGNTLYVYGEDMTPTLLRGAF

靶点信息

研究背景	The protein encoded by this gene is part of a complex termed negative elongation factor (NELF) which represses RNA polymerase II transcript elongation. This protein bears similarity to nuclear RNA-binding proteins; however, it has not been demonstrated that this protein binds RNA. The protein contains a tract of alternating basic and acidic residues, largely arginine (R) and aspartic acid (D). The gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6.
------	---

基因ID	7936
基因名	NELFE
Swiss	P18615
别名	NELFE;D6S45;NELF-E;RD;RDBP;RDP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn