

SPHK2 Rabbit pAb

货号: B17347

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	62kDa/65kDa/69kDa/80kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cytoplasm,Lysosome membrane,Membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human SPHK2 (NP_001191088.1).
序列	MNGHLEAEEQQDQRPDQELTGSWGHGPRSTLVRAKAMAPPPPLAASTPLLHGEFGSYPARGPRFALTLTSQLHIQRLR PKPEARPRGGLVPLAEVSGCCTLRSRSPSDSAAYFCIYTYPRGRRGARRRATRTFRADGAATYEENRAEAQRWATALTCLL RGLP

靶点信息

研究背景	This gene encodes one of two sphingosine kinase isozymes that catalyze the phosphorylation of sphingosine into sphingosine 1-phosphate. Sphingosine 1-phosphate mediates many cellular processes including migration, proliferation and apoptosis, and also plays a role in several types of cancer by promoting angiogenesis and tumorigenesis. The encoded protein may play a role in breast cancer proliferation and chemoresistance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	---

基因ID	56848
基因名	SPHK2
Swiss	Q9NRA0
别名	SPHK2;SK 2;SK-2;SPK 2;SPK-2;SK2;SPK2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn