

WLS Rabbit pAb

货号: B17392

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	62kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse kidney,Mouse lung,Mouse brain,Rat kidney
细胞定位	cytoplasmic vesicle,cytosol,early endosome,endoplasmic reticulum,extracellular exosome,Golgi apparatus,plasma membrane,trans-Golgi network
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 43-232 of human WLS (NP_001002292.3).
序列	YMSVKVDARKNHHKTKWFVPWGPNHCDKIRDIEEAIPREIEANDIVFSVHIPLPHMEMSPWFQFMLFILQLDIAFKLNNQIRENAEVSMDVSLAYRDDAFAEWTEMAHERVPRKLKCTFTSPKTPEHEGRYYECDVLPFMEIGSVAHKFYLLNIRLPVNEKKKINVGIGEIKDIRLVGIHQNGGFTKVW

靶点信息

研究背景	Regulates Wnt proteins sorting and secretion in a feedback regulatory mechanism. This reciprocal interaction plays a key role in the regulation of expression, subcellular location, binding and organelle-specific association of Wnt proteins. Plays also an important role in establishment of the anterior-posterior body axis formation during development (By similarity).
------	--

基因ID	79971
基因名	WLS
Swiss	Q5T9L3
别名	EVI;MRP;GPR177;mig-14;C1orf139;WLS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn