

CD215/IL-15R alpha Rabbit pAb

货号: B17398

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	15-28kDa
实测分子量	60KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1
细胞定位	Cytoplasmic vesicle membrane,Endoplasmic reticulum membrane,Golgi apparatus membrane,Membrane ,Nucleus membrane,Secreted,Single-pass type I membrane protein,Single-pass type I membrane protein, extracellular space
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-205 of human CD215 /IL-15R alpha (NP_002180.1).
序列	ITCPPPMSVEHADIWVKSYSLSRERYICNSGFKRKAGTSSLTECVLNKATNVAHWTTPSLKCIRDPALVHQRPAAPPSTVTT AGVTPQPESLSPSGKEPAASSPSSNNTAATTAAIVPGSQLMPSPKSPSTGTTEISSHESSHGTPSQTTAKNWELTASASHQP PGVYPQGHSDTT

靶点信息

研究背景	This gene encodes a cytokine receptor that specifically binds interleukin 15 (IL15) with high affinity. The receptors of IL15 and IL2 share two subunits, IL2R beta and IL2R gamma. This forms the basis of many overlapping biological activities of IL15 and IL2. The protein encoded by this gene is structurally related to IL2R alpha, an additional IL2-specific alpha subunit necessary for high affinity IL2 binding. Unlike IL2RA, IL15RA is capable of binding IL15 with high affinity independent of other subunits, which suggests distinct roles between IL15 and IL2. This receptor is reported to enhance cell proliferation and expression of apoptosis inhibitor BCL2L1/BCL2-XL and BCL2. Multiple alternatively spliced transcript variants of this gene have been reported.
基因ID	3601
基因名	IL15RA
Swiss	Q13261
别名	IL15RA;CD215

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn