

PROKR2 Rabbit pAb

货号: B17411

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	43kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,SKOV3,293T,Mouse thymus
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human PROKR2 (NP_658986.1).
序列	MAAQNGNTSFTPNFNPPQDHASSLSFNFSYGDYDLPMDDEDMDTKTRTFFAAKIVIGIAL

靶点信息

研究背景	Prokineticins are secreted proteins that can promote angiogenesis and induce strong gastrointestinal smooth muscle contraction. The protein encoded by this gene is an integral membrane protein and G protein-coupled receptor for prokineticins. The encoded protein is similar in sequence to GPR73, another G protein-coupled receptor for prokineticins.
基因ID	128674

基因名	PROKR2
Swiss	Q8NFJ6
别名	PROKR2;GPR73L1;GPR73b;GPRg2;HH3;KAL3;PKR2;dj680N4.3

产品验证



Western blot analysis of PROKR2 expressed in Mouse thymus,U-251MG,SKOV3 using PROKR2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn