

FAF1 Rabbit pAb

货号: B17423

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	56kDa/73kDa
实测分子量	85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR3,HT-29,Mouse brain,Mouse liver,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human FAF1 (N P_008982.1).
序列	MASNMDREMILADFQACTGIENIDEAITLLEQNNWDLVAAINGVIPQENGILQSEYGGETIPGPAFNPASHPASAPTSSSSS AFRPVMPSRQIVERQPRMLDFRVEYRDRNVDVLEDTCVGEIKQILENELQIPVSKMLLKGWKTGDVEDSTVLKSLHLPK NNSLYVLTPDLPPSSSSSHAGALQESLNQNFMLIITHREVQREYNLNFSGSSTIQEVKRNVDLTIPVRHQLWEGWPTSA TDDSMCLAESGLSYPCHRLTVGRRSSPAQTREQSEE

靶点信息

研究背景	Interaction of Fas ligand (TNFSF6) with the FAS antigen (TNFRSF6) mediates programmed cell death, also called apoptosis, in a number of organ systems. The protein encoded by this gene binds to FAS antigen and can initiate apoptosis or enhance apoptosis initiated through FAS antigen. Initiation of apoptosis by the protein encoded by this gene requires a ubiquitin-like domain but not the FAS-binding domain.
------	--

基因ID	11124
基因名	FAF1
Swiss	Q9UNN5
别名	FAF1;CGI-03;HFAF1s;UBXD12;UBXN3A;hFAF1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn