

CFHR1 Rabbit pAb

货号: B17460

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	37kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa,NCI-H460,Mouse liver,Mouse lung,Mouse heart,Rat liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 151-330 of human CFHR1 (NP_002104.2).
序列	PTVQNAHILSRQMSKYPSEGERVRYECSRPYEMFGDEEVMCLNGNWTEPPQCKDSTGKCGPPPIDNGDITSFPLSVYAPASSVEYQCQONLYQLEGNKRITCRNGQWSEPPKCLHPCVISREIMENYNIALRWTAKQKLYLRTGESAEFVCKRGYRLSSRSH TLRTTCWDGKLEYPTCAKR

靶点信息

研究背景	This gene encodes a secreted protein belonging to the complement factor H protein family. It binds to Ps eudomonas aeruginosa elongation factor Tuf together with plasminogen, which is proteolytically activate d. It is proposed that Tuf acts as a virulence factor by acquiring host proteins to the pathogen surface, co ntrolling complement, and facilitating tissue invasion. Mutations in this gene are associated with an incre ased risk of atypical hemolytic-uremic syndrome.
------	--

基因ID	3078
基因名	CFHR1
Swiss	Q03591
别名	CFHR1;CFHL;CFHL1;CFHL1P;CFHR1P;FHR1;H36-1;H36-2;HFL1;HFL2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn