

FES Rabbit pAb

货号: B17491

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	78kDa/85kDa/86kDa/93kDa
实测分子量	98kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse brain
细胞定位	Cell junction,Cell membrane,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle,Golgi apparatus,Peripheral membrane protein,cytoskeleton,cytosol,focal adhesion
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-420 of human FES (NP_001137256.1).
序列	WQQLQQELTKTHSQDIEKLKSQYRALARDSAQAKRKYQEASKDKDRDKAKDKYVRSWLKLFahhNRYVLGVRAAQLHH QHHHQLLLPGLLRSLQDLHEEMACILKEILQEYLEISSLVQDEVVAIHREMAAAAARIQPEAEYQGFLRQYGSAPDVPPCVTF DESLLEEGEPLPGEQLNELTVESVQHTLTSTVDELAVATEMVFRREQEMVTQLQQELRNEEENTHPRERVQLLGKRQVLQ EALQGLQVALCSQAKLQAQQELLQTKLEHLGPGEPPPVLQLQDDRHSSTSSSEQEREGR

靶点信息

研究背景	This gene encodes the human cellular counterpart of a feline sarcoma retrovirus protein with transforming capabilities. The gene product has tyrosine-specific protein kinase activity and that activity is required for maintenance of cellular transformation. Its chromosomal location has linked it to a specific translocation event identified in patients with acute promyelocytic leukemia but it is also involved in normal hematopoiesis as well as growth factor and cytokine receptor signaling. Alternative splicing results in multiple variants encoding different isoforms.
基因ID	2242
基因名	FES
Swiss	P07332
别名	FES;FPS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn