

IRF2 Rabbit pAb

货号: B17498

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	39kDa
实测分子量	45KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse lung
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-349 of human IRF2 (NP_002190.2).
序列	KKGKKPKTEKEDKVKHIKQEPVESSLGLSNGVSDLSPYAVLTSTIKNEVDSTVNIIVVGQSHLDSNIENQEIVTNPPDICQV VEVTTESDEQPVMSELYPLQISPVSSEYAESETTDSVPSDEESAEGRPHWRKRNIIEGKQYLSNMGTRGSYLLPGMASFVT SNKPDQLQVTIKESNPVPYNSSWPPFQDLPLSSSMTPASSSSRPDRETRASVIKKTSDITQARVKSC

靶点信息

研究背景	IRF2 encodes interferon regulatory factor 2, a member of the interferon regulatory transcription factor (IRF) family. IRF2 competitively inhibits the IRF1-mediated transcriptional activation of interferons alpha and beta, and presumably other genes that employ IRF1 for transcription activation. However, IRF2 also functions as a transcriptional activator of histone H4.
------	--

基因ID	3660
基因名	IRF2
Swiss	P14316
别名	IRF2;IRF-2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn