

EFNB1 Rabbit pAb

货号: B17511

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	38kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,NCI-H460,Mouse lung,Mouse kidney
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 28-237 of human EFNB1 (NP_004420.1).
序列	LAKNLEPVSWSLNPKFLSGKGLVIYPKIGDKLDIICPRAEAGRPYEYKLYLVLRPEQAAACSTVLDPNVLVTCNRPEQEIRFTIKFQEFSPNYMGLFCKKHDYYITSTSNGSLEGLNREGGVCRTRTMKIIMKVGQDPNAVTPQLTTSRPSKEADNTVKMATQAPGSRGSLGSDGKHETVNQEEKSGPGASGGSSGDPDGGFFNSK

靶点信息

研究背景	The protein encoded by this gene is a type I membrane protein and a ligand of Eph-related receptor tyrosine kinases. It may play a role in cell adhesion and function in the development or maintenance of the nervous system.
基因ID	1947

基因名	EFNB1
Swiss	P98172
别名	EFNB1;CFND;CFNS;EFB1;EFL3;EPLG2;Elk-L;LERK2;ephrin-B1

产品验证



Western blot analysis of EFNB1 expressed in Mouse lung,HT-29,NCI-H460 using EFNB1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn