

ATP2C1 Rabbit pAb

货号: B17513

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	96- 107kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-431,Rat brain
细胞定位	Golgi apparatus membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 400-660 of human ATP2 C1 (NP_001186114.1).
序列	AVSRIVEAGCVCNDAVIRNNTLMGKPTEGALIALAMKMGLDGLQQDYIRKAEYPFSSEQKWMVAVKCVHRTQQDRPEICFM KGAYEQVIKYCTTYQSKGQTLTLTQQQRDVYQQEKARMGSAGLRVLALASGPQLTFLGLVGIIDPPRTGVKEAVTTLIA SGVSIKMITGDSQETAVAIASRLGLYSKTSQSVSGEEIDAMDVQQLSQIVPKVAVFYRASPRHKMKIISLQKNGSVVAMT GDGVNDAVALKAADIGVA

靶点信息

研究背景	The protein encoded by this gene belongs to the family of P-type cation transport ATPases. This magnesium-dependent enzyme catalyzes the hydrolysis of ATP coupled with the transport of calcium ions. Defects in this gene cause Hailey-Hailey disease, an autosomal dominant disorder. Alternatively spliced transcript variants encoding different isoforms have been identified.
------	--

基因ID	27032
基因名	ATP2C1
Swiss	P98194
别名	ATP2C1;ATP2C1A;BCPM;HHD;PMR1;SPCA1;hSPCA1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn