

MTA3 Rabbit pAb

货号: B17539

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	58kDa/67kDa
实测分子量	67kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 40-330 of human MTA3 (NP_065795.1).
序列	AKVVCFYRRRDISNTLIMLADKHAKIEEESETTVEADLTDKQKHQLKHRELFSLRQYESLPATHIRGKCSVALLNETESVLS YLDKEDTFFYSLVYDPSLKTLLADKGEIRVGPRYQADIPEMLLEGESDEREQSKLEVKVWDPNSPLTDRQIDQFLVVARAV GTFARALDCSSSVRQPSLHMSAAAAASRDITLFHAMDTLYRHSYDLSSAISVLVPLGGPVLCRDEMEEWSASEASLFEEALE KYGKDFNDIRQDFLPWKSLTSIIEYYMWKTTDRYVQQKRLKAAEA

靶点信息

研究背景	Plays a role in maintenance of the normal epithelial architecture through the repression of SNAI1 transcrip tion in a histone deacetylase-dependent manner, and thus the regulation of E-cadherin levels. Contribut es to transcriptional repression by BCL6.
基因ID	57504

基因名	MTA3
Swiss	Q9BTC8
别名	MTA3

产品验证



Western blot analysis of MTA3 expressed in A-549 using MTA3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn