

PRMT4/CARM1 Rabbit pAb

货号: B17548

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200
理论分子量	41kDa/63kDa/65kDa
实测分子量	63KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 486-603 of human PRMT4/CARM1 (NP_954592.1).
序列	GSHYTSPENMWNTGSTYNLSSGMAVAGMPTAYDLSSVIASGSSVGHNLIPLANTGIVNHTHSRMGSIMSTGIVQGSSGAQGSGGGGSAHYAVNSQFTMGGPAISMASPMSIPTNT

靶点信息

研究背景	This gene belongs to the protein arginine methyltransferase (PRMT) family. The encoded enzyme catalyzes the methylation of guanidino nitrogens of arginyl residues of proteins. The enzyme acts specifically on histones and other chromatin-associated proteins and is involved in regulation of gene expression. The enzyme may act in association with other proteins or within multi-protein complexes and may play a role in cell type-specific functions and cell lineage specification. A related pseudogene is located on chromosome 9.
------	---

基因ID	10498
基因名	CARM1
Swiss	Q86X55
别名	CARM1;PRMT4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn