

LETMD1 Rabbit pAb

货号: B17565

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	11kDa/17kDa/23kDa/35kDa/41kDa/43kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse liver
细胞定位	Mitochondrion outer membrane,Single-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 161-360 of human LETMD1 (NP_056231.3).
序列	FPRQLLIRHFWTPKQQTDFLDIYHAFRKQSHPEIISYLEKVIPLISDAGLRWRLTDLCTKIQRGTHPAIHDILALRECFSNHPLGMNQLQALHVKALSRAMLLTSYLPPLLRHRLKTHTTVIHQDKALAKLGIGQLTAQEVKSACYLRGLNSTHIGEDRCRTWLGEWLQISCSLKEAELSLLLHNVVLLSTNYLGTRR

靶点信息

研究背景	This gene encodes a mitochondrial outer membrane protein. It has a potential role in tumorigenesis, which may result from negative regulation of the p53 tumor suppressor gene. Alternatively spliced transcript variants have been noted for this gene.
------	--

基因ID	25875
基因名	LETMD1
Swiss	Q6P1Q0
别名	LETMD1;1110019O13Rik;HCCR;HCCR-1;HCCR-2;HCCR1;HCCR2;HCRR-2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn