

# GNA15 Rabbit pAb

货号: B17580

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Mouse,Rat                                                                                            |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000                                                                            |
| 理论分子量 | 43kDa                                                                                                |
| 实测分子量 | 43KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse spleen,Mouse lung,Rat spleen                                                                   |
| 细胞定位  | plasma membrane,synapse                                                                              |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                                                                      |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 112-374 of human GNA 15 (NP_002059.3).                                                                                                                                                                 |
| 序列   | HASLVMSQDPYKVTTFEKRYAAAMQWLWRDAGIRACYERRREFHLLDSAVYYLSHLERITEEGYVPTAQDVLRSRMPPTTGI<br>NEYCFSVQKTNLRIVDVGGQKSERKKWIHCFENVIALIYLASLSEYDQCLEENNQENRMKESLALFGTILELPWFKSTSVILF<br>LNKTDILEEKIPTSHLATYFPSFQGPQDAEAAKRFILDMYTRMYTGCVDGPEGSKKGARSRRLLFSHYTCATDTQNIRKVFK<br>DVRDSVLARYLDEINLL |

靶点信息

|      |                                                                                                                                        |
|------|----------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Guanine nucleotide-binding proteins (G proteins are involved as modulators or transducers in various tran smembrane signaling systems. |
| 基因ID | 2769                                                                                                                                   |

|       |             |
|-------|-------------|
| 基因名   | GNA15       |
| Swiss | P30679      |
| 别名    | GNA15;GNA16 |

产品验证



Western blot analysis of GNA15 expressed in Mouse spleen,Mouse lung,Rat spleen using GNA15 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)