

DGKZ Rabbit pAb

货号: B17650

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	101kDa/103kDa/104kDa/106kDa/124kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa,293T,Mouse brain
细胞定位	Cell membrane,Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 808-1117 of human DGKZ (NP_001099010.1).
序列	ATMVQKAKRRSAAPLHSDQQPVPEQLRIQVSRVSMHDYEALHYDKEQLKEASVPLGTVVVPGDSLELCRAHIERLQQEPDGAGAKSPTCQKLSPKWCFLDATTASRFYRIDRAQEHLNYVTEIAQDEIYILDPELLGASARPDLPPTSPLPTSPCSPTPRS LQGDAAPPQGEELIEAAKRND FCKLQELHRAGGDL MHRDEQSR TLLHHAVSTGSKDVVRYLLDHAPPEILD AVEENGETC LHQAAALGQRTICHYIVEAGASLMKTDQQGDTPRQRAEKAQDTELAAYLENRQHYQMIQREDQETAV

靶点信息

研究背景	The protein encoded by this gene belongs to the eukaryotic diacylglycerol kinase family. It may attenuate protein kinase C activity by regulating diacylglycerol levels in intracellular signaling cascade and signal transduction. Alternative splicing occurs at this locus and multiple transcript variants encoding distinct isoforms have been identified.
------	---

基因ID	8525
基因名	DGKZ
Swiss	Q13574
别名	DGKZ;DAGK5;DAGK6;DGK-ZETA;hDGKzeta

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn