

BCLAF1 Rabbit pAb

货号: B17663

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	85kDa/100kDa/105kDa/106kDa
实测分子量	150KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y
细胞定位	Cytoplasm,Nucleus,Nucleus speckle,nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 600 to the C-terminus of human BCL AF1 (NP_001070909.1).
序列	REEEFKTHHEMKEYSGFAGVSRPRGTFFRIRGRGRARGVFAGTNTGPNNSNTTFQKRPKEEEWDP EYTPKSKKYFLHDD RDDGVDYWAKRGRGRGTFQRGRGRFNFKKSGSSPKWTHDKYQGDGIVEDEEETMENNEEKKDRRKEEKE

靶点信息

研究背景	This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of p roteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of B CL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zo ne near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding differen t isoforms have been found for this gene.
基因ID	9774

基因名	BCLAF1
Swiss	Q9NYF8
别名	BCLAF1;BTF;bK211L9.1

产品验证



Western blot analysis of BCLAF1 expressed in SH-SY5Y using BCLAF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn