

# SOX21 Rabbit pAb

货号: B17669

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	28kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 167-276 of human SOX 21 (NP_009015.1).
序列	GSPYSLDLGSKMAEISSSSSGLPYASSLGYPYTAGAGAFHGAAAAAAAAAAAAAGGHTSHSPGPNPGYMIPCNCSAWPS PGLQPPLAYILLPGMGKPQLDPYPAAYAAAL

## 靶点信息

研究背景	SRY-related HMG-box (SOX) genes encode a family of DNA-binding proteins containing a 79-amino acid HMG (high mobility group) domain that shares at least 50% sequence identity with the DNA-binding HMG box of the SRY protein (MIM 480000). SOX proteins are divided into 6 subgroups based on sequence similarity within and outside of the HMG domain. For additional background information on SOX genes, see SOX 1 (MIM 602148).
基因ID	11166

基因名	SOX21
Swiss	Q9Y651
别名	SOX21;SOX25

产品验证



Western blot analysis of SOX21 expressed in Mouse kidney using SOX21 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)