

# LZTS1 Rabbit pAb

货号: B17671

产品信息

|       |                                                                                                                               |
|-------|-------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                               |
| 宿主    | Rabbit                                                                                                                        |
| 克隆性   | Polyclonal                                                                                                                    |
| 预测反应  |                                                                                                                               |
| 应用    | WB                                                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:1000 - 1:5000                                                                                                    |
| 理论分子量 | 8kDa/22kDa/55kDa/59kDa/62kDa/63kDa/66kDa                                                                                      |
| 实测分子量 | 67kDa                                                                                                                         |
| 形式    | Liquid                                                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.                          |
| 偶联物   | Unconjugated                                                                                                                  |
| 阳性对照  | Jurkat,Mouse brain,Rat brain                                                                                                  |
| 细胞定位  | Cell junction,Cell membrane,Cell projection,Cytoplasm,dendritic spine,postsynaptic cell membrane,postsynaptic density,synapse |
| 纯化    | Affinity purification                                                                                                         |


抗原信息

|      |                                                                                                                                                                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 337-596 of human LZTS 1 (NP_066300.1).                                                                                                                                                 |
| 序列   | QLQQEKRQLRQELESLMKEQDLLETKLRSYEREKTSFGPALEETQWEVCQKSGEISLLKQQLKESQTEVNAKASEILGLKAQLKDTRGKLEGLELRTQDLEGALRTKGLELEVCELENLQRKKNEAELLREKVNILLEQELQELRAQAALARDMGPPTFPEDVPALQRELERLRAELREERQGHQMSSGFQHERLVWKEEKEKVIQYQKQLQQSYVAMYQRNQRLEKALQQLARGDSAGEPLEVDLEGADIPYEDIATEI |

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene encodes a tumor suppressor protein that is ubiquitously expressed in normal tissues. In uveal melanomas, expression of this protein is silenced in rapidly metastasizing and metastatic tumor cells but has normal expression in slowly metastasizing or nonmetastasizing tumor cells. This protein may have a role in cell-cycle control by interacting with the Cdk1/cyclinB1 complex. This gene is located on chromosomal region 8p22. Loss of heterozygosity (LOH) in the 8p arm is a common characteristic of many types of cancer. |
| 基因ID  | 11178                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 基因名   | LZTS1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Swiss | Q9Y250                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| 别名    | LZTS1;F37;FEZ1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)