

ZNF182 Rabbit pAb

货号: B17676

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	71kDa/73kDa
实测分子量	74kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-210 of human ZNF182 (NP_001171570.1).
序列	QVDEQIERQHQDDQDKCLLMQVGFSDDKTIITKSARDCHEFGNHLHLSTNLVASIQRPDKHESFGNNMVDNLDLFSRSSA ENKYDNGCAKLFFHTEYEKTNPGMKPYGYKE

靶点信息

研究背景	Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This gene encodes a zinc finger protein, and belongs to the krueppel C2H2-type zinc-finger protein family. It may be involved in transcriptional regulation. Multiple alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	7569

基因名	ZNF182
Swiss	P17025
别名	ZNF182;HHZ150;KOX14;ZNF21;Zfp182

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn